



Dandle PAL™

Adaptable Boundary with Innovative Dual Weighted Pouches



Dual pouch design provides more height, support and stability where it is needed most

Strap helps to stabilize the smaller pouch on top of the larger pouch by pulling it under the 'arms'



'Arms' are extra long to provide support and flexion through the trunk. Longer 'arms' also give ample room to move the micro discs around to individualize support for the infant

The Dandle PAL is made up of two sections or pouches that are filled with micro discs that feel silky to the touch. The weighted pouches provide flexible boundaries that can be easily molded around the infant to provide support.

- o NEW! Two Dandle PAL options
 - Wipeable with disposable or reusable cloth covers
 - Washable
- o Dual pouch design offers flexibility
 - Larger pouch fills the 'arms' that provide support and flexion
 - Smaller pouch gives height or extra support where needed
 - Micro discs can be moved around within each pouch to individualize care
- o Serves a variety of functions
 - Acts like an extra set of hands for comforting baby during patient care
 - Holds medical devices in place to reduce the risk of malposition, unintended removal or pressure injury
 - Supports upper and lower extremities simultaneously when used in pairs
- o May be used in conjunction with other positioning aids in the DandleLION Neurodevelopmental Support System

Dandle PAL



Two Dandle PALs are used to provide containment and boundaries entirely around the baby. Notice the 'arms' are long enough to help round the shoulders and keep the hips in flexion.



The Dandle PAL keeps the head in midline while providing support to the trunk and rounding the shoulders. Pouches are stacked one on top of the other and secured with the strap, ensuring a firm, stable surface to support the CPAP tubing at the desired height.



The PAL around the baby's head supports the CPAP equipment and the double pouches provide firm support behind the head to keep the baby from arching.



The smaller pouch is folded on top of the large pouch to provide firmness and height, which helps keep baby snugly in place without sliding down the inclined mattress.

Ordering Information

D 11505 B10	Dandle PAL- Washable	10/case
D 20505 B10	Dandle PAL- Wipeable, with white disposable cover	10/case
D 20515 C B25	Reusable printed cotton covers for DandlePAL- Wipeable	25/case
D 20525 C B25	White disposable covers for DandlePAL- Wipeable	25/case