



Dandle ROO™ Lite

Position and Support Babies During Phototherapy

Modular ROO accessories can be added to meet your needs

- Add adjustable brim to help shield eyes from direct light
- Add multi-purpose roll as a prone roll or supine support to bring hips into a posterior tilt
- Add gel pillow and cotton cover for a soft supportive surface
- All Dandle ROO products can be used in a variety of ways to meet the baby's needs



Cost effective, Dandle ROO Lite can be used during phototherapy

- Made from soft, stretchy cotton, not synthetic non-woven material
- Lightweight material allows at least 70% of phototherapy light to penetrate through a single layer
- May be hand laundered or just thrown away if laundry facilities are not available!
- Adjustable micro-fasteners provide secure support
- Upper and lower extremities can be adjusted independently of each other and may be individualized according to the baby's medical or developmental needs
- Headroll may be fastened out of the way to accommodate and/or support medical equipment
- Think creatively!! Use for admissions, "messy" procedures, eye exams, hearing screens, during transport, for agitated babies during phototherapy and any time a baby would benefit from containment or calming!
- Comes in four sizes: XS for babies <1000 g (D 11031 B12, 12/case); S for babies 1000 – 1800 g (D 11032 B12, 12/case); M for babies 1800 – 2500 g (D 11033 B12, 12/case); and L for babies 2500 – 3600 g (D 11034 B12, 12/case).

