



# Dandle•LION Polly PAL™

Versatile Positioning Support Filled with Micro Discs

## About the Polly PAL:

- Single bag design
- Filled with silky micro discs
- Helps achieve and maintain specific positions of baby or equipment
- May be placed around head, neck or hips
- May be partially placed on top of the baby for comfort  
**Caution: Do not place the entire Polly PAL on top of the baby**
- Washable according to CDC guidelines
- Reusable between patients

Polly PAL, little sister to the Dandle PAL, is a small, flexible boundary that supports positioning of NICU babies.

Each Polly PAL is filled with tiny plastic micro discs, giving it a silky feel. The Polly PAL can be used in a myriad of ways, including:

- Around the head to encourage midline positioning
- Supporting respiratory equipment or other medical tubing
- Behind the back when baby is side-lying
- As a “second set of hands” to provide comfort during care times
- Anywhere the baby needs a little more comfort or support



## Ordering Information

D 21755 B10 | Polly PAL Infant Positioning Aid | Qty: 10/case

Made in USA

Bimeco Group, INC, P: 770.631.7761, F: 770.631.7793, [customerservice@bimecogroup.com](mailto:customerservice@bimecogroup.com)

**Dandle•LION™**  
MEDICAL

Rev1:0621